## DISPOSABLE CRANIAL PERFORATORS INSTRUCTIONS FOR USE

English  
Czech/Čeština  
Danish/Dansk  
Dutch/Nederlands  
French/Français  
German/Deutsch  
Greek/Ελληνικά  
Hungarian/Magyar  
Italian/Italiano  
Lithuanian/ Lietuvių  
Norwegian/Norsk  
Polish/Polski  
Portuguese/Português  
Spanish/Español  
Swedish/Svenska  
Japanese/日本語  
Chinese/中文

**DISPOSABLE CRANIAL PERFORATOR**

**For Single Patient Use**

**FOR THE FOLLOWING MODELS:**

<b>Model 200-331</b>	<b>11/7mm</b>	<b>DG-II Mini</b>
<b>Model 200-251</b>	<b>14/11mm</b>	<b>DGR-II</b>
<b>Model 200-253</b>	<b>11/7mm</b>	<b>DGR-II Mini</b>
<b>Model 200-255</b>	<b>13/9mm</b>	<b>DGR-II</b>
<b>Model 200-241</b>	<b>14/11mm</b>	<b>DGR-I</b>
<b>Model 200-243</b>	<b>11/7mm</b>	<b>DGR-I Mini</b>
<b>Model 200-245</b>	<b>13/9mm</b>	<b>DGR-I</b>
<b>Model 200-500</b>	<b>14/9mm</b>	<b>SDR</b>
<b>Model 200-271</b>	<b>14/11mm</b>	<b>DGR-0</b>
<b>Model 200-275</b>	<b>13/9mm</b>	<b>DGR-0</b>
<b>Model 200-283</b>	<b>8/5mm</b>	<b>DGR-0S</b>
<b>Model 210-221</b>	<b>14/11mm</b>	

**IMPORTANT  
INFORMATION  
Read before using  
FAILURE TO READ  
AND FOLLOW THESE  
INSTRUCTIONS MAY  
RESULT IN SERIOUS  
PATIENT INJURY,  
DEATH OR PRODUCT  
MALFUNCTION.**

**Please note:** These instructions include use and cautions of the DGR-I, DGR-II, DGR-0, DGR-0S, SDR, DG-II and 210-221 Disposable Perforators. Use caution when selecting appropriate model for use in specific skull areas and skull conditions.

**CAUTION:** ACRA-CUT cranial perforators are very sharp, and therefore cut and respond differently than other brand perforators. To ensure safe and reliable performance, carefully read and follow instructions. Failure to follow instructions for use may result in serious patient injury, death or product malfunction.

**210-221 Perforators**

All Model 210-221 disposable perforators are designed to automatically release and stop upon penetration of bone that is **at least 2.5mm thick**. Caution should be taken to avoid use in the temporal, suboccipital areas and in children and/or infants.

**DG-II Perforators**

All Model DG-II disposable perforators are for use on thin skull/skull areas such as pediatric, temporal, and suboccipital areas. They are designed to automatically release and stop upon penetration of bone **as thin as 1mm**. Caution should be taken to avoid excessive force when drilling through thin skulls to prevent pushing the entire unit through the skull.

**DGR-II Perforators**

All Model DGR-II disposable perforators are for use on thin skull/skull areas such as pediatric, temporal, and suboccipital areas. They are designed to automatically release and stop upon penetration of bone **as thin as 1mm**. Caution should be taken to avoid excessive force when drilling

through thin skulls to prevent pushing the entire unit through the skull.

#### **SDR Perforators**

All Model SDR disposable perforators are designed to automatically release and stop upon penetration of bone that is **at least 1.5mm thick**. Caution should be taken to avoid use on children and/or infants.

#### **DGR-0S Perforators**

All Model DGR-0S disposable perforators are designed to automatically release and stop upon penetration of bone that is **at least 2mm thick**. Caution should be taken to avoid use on children and/or infants.

#### **DGR-I and DGR-0 Perforators**

All Model DGR-I and DGR-0 disposable perforators are designed to automatically release and stop upon penetration of bone that is **at least 3mm thick**. Caution should be taken to avoid use in the temporal, suboccipital areas and in children and/or infants.

#### **ALL MODELS**

##### **1.0 Single Patient Use:**

- 1.1** This perforator has been shipped sterile and is for single patient use only.

- 1.2** Resterilization will cause serious damage to the unit and therefore to the patient.

##### **2.0 Testing Procedure:**

Unit has been presterilized. Handle accordingly. Perform the following test before each perforation. Protect gloved fingers with gauze while performing tests. If either test does not perform as indicated, take special action to prevent further use of the perforator and immediately return to ACRA-CUT, Inc. To test the mechanical integrity of the assembled unit:

- 2.1** Hold the drive shank securely while pressing on the fluted components (inner and outer drill). A smooth, positive spring action should be felt.
- 2.2** While holding the drive shank, rotate the fluted portion (inner and outer drill). These parts should turn smoothly.

##### **3.0 General Operating Procedures:**

- 3.1** The drive shank fits a standard Hudson coupler. Proper care should be taken to insure the drive shank is fully seated in the drive unit.
- 3.2** Visually examine the perforator and its packaging for any obvious defects, damage, or distortion. Read and inspect the label on the perforator. Do not

use perforator if brown stripe is visible on perforator label as this indicates the perforator has been resterilized by ethylene oxide, contrary to the labeling.

- 3.3 Power unit must be used only in the CLOCKWISE direction, or forward mode.
- 3.4 Perforator can be used at low speeds and in hand-operated drilling devices.
- 3.5 Recommended speed for power drilling is up to 800 RPM. Power unit (standard 1,000 RPM) should be first set at 60 PSI and increased as needed to no more than 80 PSI.
- 3.6 To start perforation, place the nonskid tip against the skull, compress and rotate drill by hand until drive mechanism is seated in the driver. Then drill at chosen speed.
- 3.7 Perforator will stop when penetration is made or when surgeon releases pressure on the unit.
- 3.8 In the event of any unsatisfactory results, do not continue to use perforator.
- 3.9 Repeat steps 2.0, 2.1 and 2.2 of Testing Procedure and steps 3.2 through 3.9 of the General Operating Procedure for each perforation.

### 3.10 Shelf Formation:

**3.10.1** The shelf formed by DGR-I, DGR-0, DGR-0S, and SDR perforators will be well defined. The DGR series of perforators also form a ridge of bone at the bottom of the through hole (See Diagram - Figures 1, 3 and 4).

**3.10.2** The shelf formed by DGR-II perforators slants more than the shelf formed by DGR-I, DGR-0, and SDR perforators. This feature allows for safe use on skull/skull areas as thin as 1mm. The DGR series of perforators also form a ridge of bone at the bottom of the through hole. (See Diagram - Figure 2)

**3.10.3** The shelf formed by the Model 210-221 disposable perforator will be well defined. (See Diagram - Figure 5)

**3.10.4** The shelf formed by DG-II perforators slants more than the shelf formed by the Model 210-221 perforators. This feature allows for safe use on skull/skull areas as thin as 1mm. (See Diagram - Figure 6)

**FOR DG-II, DGR-I, DGR-II,  
DGR-0, DGR-0S AND SDR****MODELS ONLY**

**3.11** The Dura-Guard “R” Series tip forms a pad of bone at the bottom of the burr hole. This pad protects the dura from the cutting edges and pushes the dura away from the skull. After irrigation and suctioning, the pad is easily removed with a pick or elevator type instrument such as the ACRA-CUT Model 800-020 Dura Separator. If it is suspected that the dura is adhering to the bone, EXTRA CAUTION must be taken or dura may be cut. (See Diagrams - Figure 7 and Figure 8)

**WARNINGS**

- 4.0** Failure to heed warnings or to follow procedures as described could result in serious patient injury, death or product malfunction.
- 4.1** The single patient use perforator must never be used on more than one patient. Multiple use could cause product malfunction, therefore serious patient injury or death.
- 4.2** Surgeon should always be aware of the possibility that the perforator may not disengage.
- 4.3** Excessive hand pressure must be avoided to preclude penetration of the dura and/or

brain. Allow the perforator to do the cutting.

- 4.4** Reduce pushing force and observe caution as the perforator nears the point of penetration as nicking, cutting, or tearing the dura can occur.
- 4.5** Hold power unit/perforator perpendicular (90°) to the skull at point of penetration and do not “rock” to speed cutting action. Failure to observe could cause product malfunction, patient injury or dural laceration.
- 4.6** Extreme caution must be taken when perforating in areas of the skull that have 1mm or greater variations in bone thickness. Drill may nick or cut dura or brain (similar effect as not being 90°).
- 4.7** If conditions such as adherent dura, high intracranial pressure or other underlying abnormalities are present in the area of the penetration, the perforator may cut the dura or brain.
- 4.8** Speed for power drilling should not exceed 800 RPM. See Paragraph 3.5 for procedure.
- 4.9** Although the perforators are designed to release safely and reliably, (DGR-I/DGR-0 units in skull/skull areas at least 3mm thick; DGR-0S in skull/skull areas at least 2mm thick; DGR-II/DG-II units in skull/

skull areas as thin as 1mm; 210-221 unit in skull areas at least 2.5mm thick and SDR units in skull areas at least 1.5mm thick) caution must be observed when perforating thin skull areas, such as temporal bone, infants, children, elderly or diseased bone, since skull consistencies and thicknesses can vary and the dura or brain could be cut.

- 4.10** Do not use perforator in or near previously drilled areas.
- 4.11** Do not use perforator if package or label is damaged or missing.
- 4.12** Dispose of perforator after use according to established hospital procedure for contaminated items.



**CAUTION:** Federal (USA) law restricts this device to use by or at the order of a physician.

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ACRA-CUT Inc.
MT Promedt Consulting GmbH