

# XPRESS 120

ACRA  
CUT®

## Disposable All-in-One Battery Powered Drill with ACRA-CUT Perforator

- Drill available with various size perforators:  
14/11mm, 14/9mm, 11/7mm, 8/5mm
- Portable/Ready for use
- Trauma, E.R., O.R., bedside, etc.
- Sterile/Single patient use
- Single or multiple burr holes

989 Main St., Acton, MA 01720

[www.acracut.com](http://www.acracut.com)

Call: 800/227-2288

In MA: 978/263-0250

Fax: 978/263-4102

**For Convenience, Safety and Reliability,  
there's only ONE choice - "THE ACRA-CUT XPRESS 120"**

**A DISPOSABLE BATTERY POWERED DRILL AND ACRA-CUT PERFORATOR**

Drill single or multiple burr holes with the safety, confidence and reliability of an ACRA-CUT perforator with the added convenience of a disposable battery powered drill.

XPRESS 120 is ideal for a variety of neurosurgical applications:

- Shunt Placements
- Ventriculostomies
- ICP Catheter Placements
- Tumor Biopsy
- Stereotactic Procedures
- Emergency Backup
- Trauma
- E.R. / O.R.
- Bedside

**Additional Benefits of XPRESS 120  
Disposable Battery Powered Drill  
and Perforator**

- Eliminates sterilization inconvenience
- No set up time required, ready to drill
- Disposable single patient use
- Perforators available in 14/11mm, 14/9mm, 11/7mm and 8/5mm sizes



Includes extra battery and charger



120-241 DGR-I  
14/11mm  
For skulls at least  
3mm thick



120-271 DGR-O  
14/11mm  
For skulls at least  
3mm thick, uneven bone



120-500 SDR  
14/9mm  
For skulls as thin as  
1.5mm thick



120-253 DGR-II  
11/7mm  
For skulls as thin as  
1mm thick



120-283 DGR-OS  
8/5mm Mini  
For skulls at least  
2mm thick

Available in the following sizes:

CATALOG #	DESCRIPTION
120-241	Disposable Battery Powered Drill with 14/11mm DGR-I Perforator
120-253	Disposable Battery Powered Drill with 11/7mm DGR-II Perforator
120-271	Disposable Battery Powered Drill with 14/11mm DGR-O Perforator
120-283	Disposable Battery Powered Drill with 8/5mm DGR-OS Perforator
120-500	Disposable Battery Powered Drill with 14/9mm SDR Perforator

**Features of all ACRA-CUT Disposable  
Cranial Perforators**

- ACRA-CUT's cammed lug release mechanism is the safest and most reliable ever developed.
- Inner drill forms a bone pad to help protect the dura.
- Rigorous quality standards, including pre-test drilling seven holes with every disposable perforator make ACRA-CUT perforators the safest in the world.
- Sharp edges cut quickly and smoothly so less hand pressure is required.
- Will not stop drilling before the burr hole is completed. Can be stopped and restarted as the surgeon requires.
- Non-skid tip prevents the perforator from skidding on the skull during start up.
- Inner drill moves forward only 1mm to automatically release.



Contact ACRA-CUT:

- To arrange a product trial.
- To place an order.
- For additional information on our complete range of cranial perforators, cutting tools and other instruments.

989 MAIN STREET, ACTON, MASSACHUSETTS 01720

or visit [www.acracut.com](http://www.acracut.com)

**Toll free 1-800-ACRA-CUT (800-227-2288) In Massachusetts (978) 263-0250 FAX (978) 263-4102**

ACRA-CUT products are manufactured under one or more of the following patents: U.S. Patent Nos. 4,600,006; 4,699,550; 4,803,982; 4,884,571; 4,951,690; 5,007,911; 5,135,532; 5,658,305; 5,665,097; U.K. Patent Nos. 2,154,484; 2,181,076; 2,184,374; 2,208,362; DE Patent Nos. 3,503,098; 3,624,860; 3,890,886; JP Patent No. 2,607,658. Other U.S. and foreign patents pending.

Printed in U.S.A. © 2004 ACRA-CUT, Inc.