



DURA SEPARATORS BONE PAD REMOVERS



Standard Dura Separators



Baker Dura Separators



Bone Pad Remover

Left
Right

ACRA-CUT DURA SEPARATORS


Novel Designs That Are Easier To Use

UNIQUE BENEFITS

- Comfortable, non-slip titanium handles allow delicate control even when wet.
- Lightweight and perfectly balanced.
- Specially shaped blades reach further under the skull so separating the dura is easier and safer.
- Available in two models and four different blade widths.
- 2.5mm and 4.0mm dura separators for small diameter pediatric burr holes.

STANDARD DURA SEPARATORS


ACRA-CUT Standard Dura Separators have a gently curved blade that allows more reach under the skull than conventional dura separators.



ACRA-CUT Standard Dura Separator compared to a conventional dura separator

BAKER DURA SEPARATORS

BAKER Dura Separators' unique blade shape provides maximum reach under the skull. The blade hugs the inner table making separation of the dura easier.



BAKER Dura Separator compared to conventional dura separator

STANDARD DURA SEPARATORS

Cat. No.	Model	Blade Width	
800-710	#0	2.5 mm	Pediatric / Small burr holes
800-711	#1	4.0 mm	Pediatric / Small burr holes
800-712	#2	6.0 mm	Similar to Penfield #3
800-713	#3	8.0 mm	Similar to Penfield #3 with scalpel guide groove.
800-714 Set of 4; includes one each #0, #1, #2, #3			

BAKER DURA SEPARATORS

Cat. No.	Model	Blade Width	
800-715	#0	2.5 mm	Pediatric / Small burr holes
800-716	#1	4.0 mm	Pediatric / Small burr holes
800-717	#2	6.0 mm	Similar to Penfield #3
800-718	#3	8.0 mm	Similar to Penfield #3 with scalpel guide groove.
800-719 Set of 4; includes one each #0, #1, #2, #3			



Left hand bone pad remover

ACRA-CUT BONE PAD REMOVERS

A convenient instrument surgeons use to remove the bone pad formed by ACRA-CUT perforators quickly and easily.



Right hand bone pad remover

UNIQUE BENEFITS

- Comfortable non-slip ribbed titanium handle with a stainless steel tip.
- Lightweight, perfectly balanced and easy to use.

BONE PAD REMOVERS

Cat. No.	
800-020 R	Bone pad remover, Right hand model
800-020 L	Bone pad remover, Left hand model

Contact ACRA-CUT:

- To arrange a product trial.
- To place an order.
- For additional information on our complete range of cranial perforators, cutting tools and other instruments.

FOR MORE INFORMATION



989 MAIN STREET, ACTON, MASSACHUSETTS 01720
or visit www.acracut.com

An educational videotape, Recent Advances in Cranial Perforation, is now available. Please call for your free copy.

Toll free 1-800-ACRA-CUT (800-227-2288) In Massachusetts (978) 263-0250 FAX (978) 263-4102

ACRA-CUT products are manufactured under one or more of the following patents: U.S. Patent No. 4,600,606; U.S. Patent No. 4,699,550; U.S. Patent No. 4,803,982; U.S. Patent No. 4,884,571; U.S. Patent No. 4,951,690; U.S. Patent No. 5,007,911; U.S. Patent No. 5,135,532; U.S. Patent No. 5,658,305; U.S. Patent No. 5,665,097; U.K. Patent No. 2,154,484; U.K. Patent No. 2,181,076; U.K. Patent No. 2,184,374; U.K. Patent No. 2,208,362; DE Patent No. 3,503,098; DE Patent No. 3,624,860; DE Patent No. 3,890,886; JP Patent No. 2,607,658. Other U.S. and foreign patents pending.

Printed in U.S.A. 4015M © 2001 ACRA-CUT, Inc.